

## GETTING STARTED GUIDE

myOxyGO<sup>®</sup>



## Before You Start

### 1 Download the App

From your mobile device, go to the **App Store (Apple)** or **Google Play (Android)**, search for the **My OxyGo App** and download the software.

Apple App Store



Google Play App Store



My OxyGo App



### 2 Plug In Your OxyGo FIT™

Connect the AC power supply cord to the concentrator. Plug the power supply AC cord into an electrical outlet.

**Do not power on the FIT concentrator.**



myOxyGO<sup>®</sup>

## Pairing


### 3 Pair Your Mobile Device With Your FIT

- ▶ Go to **Settings** on your mobile device and **turn on Bluetooth**.
- ▶ Slide the button to the right to turn on Bluetooth.
- ▶ Next, **turn on Bluetooth on your concentrator** by holding down the minus button until the Bluetooth symbol appears on the concentrator's display screen.

Go to settings



Turn on Bluetooth on phone



Turn on Bluetooth on FIT




### 4 Open My OxyGo App

Enter your provider code and click **"SUBMIT"**. The provider code can be found in your confirmation email or provided by your home care provider.

*Pairing is the term used to describe the process of wirelessly connecting two Bluetooth devices for the first time. "Pairing" allows the devices to recognize each other and creates a unique lasting connection between the specific devices.*

myOxyGO<sup>®</sup>




## 5 Search for Concentrator

Click on **“Search for Concentrator”** at the bottom of your screen.




## 6 Select Serial #

When device is found, select your **Concentrator Serial Number** on your mobile device.




## 7 Terms & Conditions

Read the **Terms and Conditions** and if you choose to accept, select the **“I Accept”** button at the bottom of your screen.



## 8 Press Bell Button

Press and hold the **bell button** on your FIT concentrator to finish pairing.




## 9 Wait to Pair

Please wait — your concentrator is pairing.

## 10 Pairing Complete

**Your OxyGo FIT is now paired with your mobile device!** Now you can turn on your FIT and use it normally. The information displayed on your FIT screen will vary depending on your concentrator’s current state.

For further information, see the My OxyGo App FAQ section at: [www.oxygo.life/app](http://www.oxygo.life/app)



**myOxyGO**

*The purpose of the My OxyGo App is to provide additional convenience to the patient. However, the User Interface Panel remains the primary source of information to which the patient should refer.*

**Location Information:** When you use our mobile applications, we receive your precise location information. We may also collect the precise location of your device when the app is running in the foreground or background / when the app is closed. We use your location information to monitor the location of our equipment, locate lost or stolen equipment, perform market analysis, and locate nearby service representatives and branch offices. We also infer your more general location information (for example, your IP address may indicate your more general geographic region). We may disclose the information we collect about you (including your location information) with homecare providers involved in your care. The My OxyGo App is compatible with the following devices: iPhone 5 and later; iPad Air; iPad Air 2; iOS 9 and later; Samsung S5 and later; Nexus 5; Nexus 6; Nexus 9; Android 5 and later. Compatible device models are subject to change at any time without notice.